

# Oligonucleotides

For Life Science Research

- > Custom **oligonucleotides**
- > **NGS** oligonucleotides
- > **RNAi** oligonucleotides
- > **Real-Time qPCR** probes
- > **Highly complex** oligonucleotides
- > **Catalogue** oligonucleotides
- > **Synthesis** reagents





## A RELIABLE EXPERIENCE

Since 1985 Eurogentec has provided high-quality reagents and custom-synthesised oligonucleotides to scientists around the globe.

## CUSTOM PRODUCTS


- All types of oligonucleotides, from 2 to 225 bases
- All chemistries: DNA, RNA, LNA®, 2'O-Me, 2'O-MOE, PNA,...
- More than 300 modifications
- All synthesis scales, from µg to grams
- Wide range of Real-Time qPCR Probes
- RNAi oligonucleotides
- Custom fill & finish

## TRUSTED QUALITY

### Optimised chemistry

- Stringent quality controls
- ISO 9001 certified quality system
- ISO 13485 certified for IVD oligonucleotides


Research, Track™ and GMP grade oligonucleotides



From lab scale to large synthesis scale



More than 300 modifications with a wide range of dyes and quenchers



ISO 9001  
ISO 13485 certified quality systems



Oli&GO™ oligonucleotide e-commerce platform



## Life Science Oligonucleotides

### Custom

### Catalogue

Custom  
oligos

NGS

RNAi

qPCR  
probes

Unique™  
oligos

Catalogue  
oligos

Synthesis  
reagents

- [p5](#) *A large range for any applications*
- [p9](#) **Custom oligos**
- [p10](#) *NGS oligos*
- [p11](#) *RNAi oligos*
- [p13](#) *qPCR probes*
- [p16](#) *Spectral properties of fluorophores & quenchers*
- [p18](#) *Unique™ oligos*
- [p19](#) **Catalogue oligos**
- [p21](#) *Synthesis reagents*
- [p22](#) *Additional services*

## Annexes

- [p23](#) *Synthesis scale vs guaranteed yield*
- [p24](#) *Shipping*
- [p24](#) *Documentation*
- [p25](#) *How to store your oligo*
- [p25](#) *How to reconstitute your oligo*
- [p25](#) *How to quantify your oligo*
- [p26](#) *IUB code*
- [p27](#) *Dyes compatibility table*
- [p28](#) *How to order*
- [p28](#) *How to pay*
- [p29](#) *How to reduce my shipping fees*
- [p29](#) *Related products*
- [p30](#) *License statements*
- [p31](#) *Trademarks and labels*

## A LARGE RANGE FOR ANY APPLICATION

Whatever your application, even for those that are most demanding (NMR, X-ray crystallography, *in vivo* animal studies...), Eurogentec can provide the highest quality oligonucleotides to meet (and exceed) your expectations!

### CUSTOM OLIGOS p.9

PCR | FISH | Pyrosequencing  
| Cloning | NMR | X-Ray  
crystallography | Mutagenesis |  
SNP Analysis

### NGS OLIGOS p.10

Next-generation sequencing

### RNAi OLIGOS p.11

Gene silencing | Antisense studies

### qPCR PROBES p.13

Real-time qPCR | Patient  
management | Diagnostic assays

### UNIQUE™ OLIGOS p.18

For highly complex oligos or if  
you don't find what you need,  
please contact us at [unique@eurogentec.com](mailto:unique@eurogentec.com)

### CATALOGUE OLIGOS p.19

Cloning | Sequencing | PNA FISH |  
qPCR Calibration

### APTAMERS p.18

Selection and Production of  
custom optimised RNA and DNA



## ⇒ BACKBONES


**Bases:** Eurogentec synthesises classic DNA and RNA based oligos but other backbones like LNA<sup>®</sup>, 2'-O-Me, 2'-O-MOE and PNA can be used to match your specifications.

**LNA<sup>®</sup>** (Locked Nucleic Acid) is a bicyclic nucleic acid with a structure locked into a rigid C3'-endo position, which favours RNA A-type helix duplex geometry. This exceptional structure confers a very strong thermal stability towards complementary DNA and RNA template suitable for hybridisation assays requiring high specificity and/or reproducibility.

**2'-O-Me RNA** are partially resistant to a variety of ribo- and deoxyribonucleases. They form more stable hybrids with complementary RNA strands than equivalent DNA/RNA sequences. They are ideal for antisense probes.

**2'-O-(2-Methoxyethyl)- oligoribonucleotides or 2'-O-MOE** have an analogue chemical structure to RNA excepted that a methoxy-ethyl residue is attached at the 2'-O-position. The chemical group at this position confers to the oligo backbone a highest nuclease resistance and a better binding affinity compared to the classical RNA molecule making it a useful tool for antisense applications.

**PNA** is an artificial DNA/RNA analogue with no charged backbone. The absence of charge repulsion between PNA and the DNA/RNA complementary strand confers a higher specificity and sensitivity. PNA are also known to be resistant to the enzyme degradation and stable over a wide range of pH, temperature and salt concentrations. These properties enable PNA to be particularly useful for FISH studies and miRNA inhibition.




## Linkages

**Phosphodiester** (default) bonds connect the 3' carbon atom and the 5' carbon of the sugar.

**Phosphorothioate** bonds possess an increased resistance against nucleases due to the substitution of a non-bridging oxygen by sulphur.

**Methylphosphonate\*** bonds are non ionic nuclease resistant linkages. Methylphosphonate/RNA duplex are not recognised by RNaseH.

\* Only available with DNA bases



## MODIFICATIONS

Oligonucleotides can be modified by **direct incorporation** during the synthesis or by **post-synthesis** labelling.


### Direct incorporation

#### 3' modifications

Since automated oligonucleotide synthesis is realised from 3' to 5', these modifications are only possible if the corresponding solid support (CPG column) is available and if the modification is compatible with the chemistries used during the synthesis. Typical examples are 3'-phosphate, 3' Biotin, 3' FAM, 3' DDQ I, 3' BHQ-1®...

#### 5' and internal modifications


Many modifications can be directly introduced at the 5' end or at internal positions of the oligonucleotides using the phosphoramidites. However these modifications need to support the somewhat harsh cleavage-deprotection conditions including a strong basic pH. Typical examples are 5' Biotin, 5' Phosphate, 5' Cholesterol, 5' FAM, 8-Oxo-dA, Biotin-dT, DABCYL-dT...



### Post-synthesis incorporation

Post-synthesis modifications may influence the yield of the reaction. A lower yield may result from poly-modifications and/or strong secondary structures. Two major post-synthesis reactions are used to introduce sensitive dyes or compounds that do not exist as phosphoramidites. In the first case the label is conjugated to an amino-modified oligonucleotide (3', 5' or on a dT) using its amino-reactive version (N-hydroxysuccinimide (NHS) ester in most cases).

The second possibility (originally also used for synthesis of molecular beacons) is the addition of a maleimide-modified label to a thiol-modified oligonucleotide.



### TABLE: MODIFICATIONS CLASSIFICATION

Fluorophores	Quenchers	Phosphates	Degenerate Bases	Non-natural bases	Modifiers	Conjugates	Biotins	Spacers
AP5 <b>AlexaFluor®</b> <b>ATTO</b> <b>BODIPY®</b> <b>Cascade Blue®</b> <b>Cy®</b> <b>DragonFly™ Orange</b> <b>DY</b> FAM Fluorescein HEX <b>HiLyte™ Fluor</b> <b>JOE</b> <b>Marina Blue®</b> <b>Oregon Green® 488</b> <b>Pacific Blue™</b> <b>Rhodamine</b> <b>ROX</b> <b>TAMRA</b> TET <b>Texas Red®</b> Yakima Yellow®	BHQ® Dabcyl Deep Dark Quencher Eclipse® Dark Quencher <b>QXL™</b> <b>TAMRA</b>	3' Phosphate 5' Phosphate	Wobbles Spikes	2' Fluoro RNA 2-Amino dA 2-Aminopurine 2'O-Me 5-Me-C 2'O-Me Propyne C, U 5,6-dihydro dU 5-Br dC, dU 5-Me dC, iso dC 7-deaza dA, dG 8-Br dA, dG 8-Oxo dA, dG AP dC C5-propyne dC, dU dA, dC, ddC, dG, dT Blocked dInosine dUracil Inverted base iso dG N4-Et dC Nitroindole O6-Me dG	Amine dR Amino Modifier Propargylamine dU Thiol Modifier Thiophosphate Triphosphate  <b>Tm Modifiers</b> LNA® PNA APdC iso dG	Acridine <b>AP conjugation</b> <b>BSA conjugation</b> Carboxy dT Cholesteryl <b>Digoxigenin</b> DNP Glyceryl <b>HRP conjugation</b> <b>Peptide conjugation</b> Psoralen <b>SBP conjugation</b>	PC-Biotin Biotin-TEG	C3 9/TEG C12 18/HEG
	<b>DID YOU KNOW?</b> Access™ Dyes from Eurogentec are a simple, customisable & cost effective solution for nucleic acid detection offering the highest performance combined to IP-friendly terms and conditions. To get more information about this service, please contact our specialists at: <a href="mailto:access@eurogentec.com">access@eurogentec.com</a>	<b>KEEP IN MIND</b> <b>IUB base codes</b> <b>B</b> = C/G/T <b>D</b> = A/G/T <b>H</b> = A/C/T <b>K</b> = G/T <b>M</b> = A/C <b>N</b> = A/C/G/T <b>R</b> = A/G <b>S</b> = C/G <b>V</b> = A/C/G <b>W</b> = A/T <b>Y</b> = C/T						

**>300**  
 MODIFICATIONS  
 AVAILABLE  
 IN DIFFERENT  
 SYNTHESIS  
 SCALES.

#### >Note

- Post-synthesis modifications are highlighted in blue.
- Some modifications can be inserted after or during the synthesis and are in red.



# CHOOSE THE RECOMMENDED PURIFICATION

The aim of any purification method is to remove the by-products resulting from the removal of the protecting groups and other synthesis by-products. To know the best purification according to each modification, consult the price list available on [www.eurogentec.com](http://www.eurogentec.com). If you are not sure which purification suits your application, then please specify "Recommended Single Purification" (additional fee) and we will choose the best purification for you.

## PURIFICATION VS APPLICATIONS

- > Isothermal sequencing
- > Cycle sequencing
- > Routine PCR
- > Hybridisation
- > DNA MicroArray
- > SNP Analysis
- > AFLP
- > OLA
- > Sensitive PCR (Diagnostic)
- > PCR Primers
- > NGS
- > *in situ* Hybridisation
- > Real-Time qPCR
- > Capillary sequencing
- > miRNA, siRNA and antisense
- In vivo* studies
- > Cloning and subcloning PCR
- > Gene synthesis
- > Gel-shift assay
- > First-strand cDNA synthesis
- > Production of cloning linkers
- > Site-directed mutagenesis
- > NMR
- > X-ray crystallography

**Classical modifications (modified bases, chemical linkers...)**

**Non radioactive labelling**

**Special modifications (G-clamp...)**  
**Labelling with fluorophores and quenchers**

**SePOP desalting** increase the purity level of the deprotected and desalted oligos up to 65-70%. It uses differential precipitation to eliminate the largest part of contaminants (truncated material < 15 bases).

**RP Cartridge•Gold** consists of a reverse phase chromatography based on the difference in hydrophobicity between the full-length product and truncated sequences. It yields to 75-80% purity. It is the best compromise for most application and the absence of residues (which may occur with HPLC) makes them suitable for cell culture uses.

**HPLC** provides a degree of purity up to 85%. Reverse Phase (RP) is based on the hydrophobic interaction of the full length oligonucleotides with alkyl chains bonded on the matrix. Ion exchange (IEX) is based on the preference of the anion - exchange resins (positively charged) for the full-length oligonucleotides.


***In vivo*** oligonucleotides are the smart choice for antisense oligonucleotides or siRNA testing at a research level before entering into pre-clinical studies. The production process of *in vivo* oligonucleotides includes the following steps: HPLC purification, desalting, sterile filtration and lyophilization.

**Polyacrylamide gel (PAGE)** separates oligonucleotides varying from only 1 base and give a purity level of 85-90%. Gel band is excised under low intensity UV. Oligonucleotide is then eluted, precipitated, quantified and packaged.

**Dual HPLC** (double RP or RP+IEX) increases the purity level up to 95%.

**UltraPureGold™** relies on a proprietary synthesis and purification process combining a synthesis on polystyrene support, special amidite, optimised deprotection and a dual purification adapted on the length and the structure of the oligos. Moreover a double quality control is performed.

## PURIFICATION VS OLIGO LENGTH




\*These values are purely indicative and only valid for an unmodified oligonucleotide of 20 bases. In addition, according to your oligonucleotides (sequences, modifications...), the purity level can be analysed by various methods (analytical HPLC, CGE...).

EUROGENTEC MANUFACTURES HIGHLY PURIFIED OLIGONUCLEOTIDES UP TO 95%.

### SYNTHESIS SCALE VS GUARANTEED YIELD

Please refer to the minimum guaranteed yield table page 23 to select the right synthesis scale or contact us at: [oligocentre@eurogentec.com](mailto:oligocentre@eurogentec.com)





# Custom oligos

Eurogentec proposes a large choice of chemistries, modifications, specifications and purifications. More than 300 modifications and several purity levels are available. ■

## SPECIFICATIONS

**Length:** From 5 to 139 bases

**Synthesis scale:** 10 nmol • 40 nmol • 200 nmol • 1000 nmol  
• 2.5 µmol • 5 µmol • 10 µmol\*

**Backbone:** DNA, RNA, LNA®, 2'-O-Me RNA, 2'-O-MOE RNA, PNA and all linkages

**Modifications:** More than 300 modifications! (see p. 7)

**Purifications:** SePOP desalting, RP-Cartridge-Gold™, HPLC, *in vivo*, PAGE, Dual HPLC, UltraPureGold™

**Quality Control:** MALDI-TOF MS

**Format:** Dried (except for unmodified SePOP desalted oligonucleotides from 15 to 39 DNA bases: 100 µM H<sub>2</sub>O by default)

**Packaging:** 2 mL tube, 96-well or 384-well plates

**Documentation:** Technical data sheet

**Shipping:** At room temperature

### DID YOU KNOW?


The combination of our vast expertise in oligonucleotide to the well reputed know how of AnaSpec in complex peptide synthesis, allows us to offer you high quality **peptide-oligo conjugates**.

### SPECIFIC NEED?

Need a high production process traceability? Discover our Track oligos on [www.eurogentec.com/track-oligos.html](http://www.eurogentec.com/track-oligos.html)

>Note

\*Larger synthesis scales are available on request.



# NGS oligos

Next-Generation Sequencing (NGS) is a high-throughput technology allowing the **massive sequencing** of nucleic acids following a DNA library preparation.

After DNA fragmentation, adapters (including indexes for multiplexing) are fused to the fragments. NGS adapters require both a **high level of purity** (no n-x side products) and the **absence of cross-contamination** (confusing index sequences).

Thanks to our long history as an oligo provider Eurogentec has developed a dedicated manufacturing process for the production of high quality NGS oligos. ■

## SPECIFICATIONS

**Quality:** Low cross-contamination (<0,1%)

**Length:** from 20 to 85 bases

**Quantity:** 10 nmol minimum delivered\*

**Purification:** HPLC or cartridge

**QC:** 100% QC checked by MalDI-TOF MS

**5' Modifications:** 5' Phosphate / 5' Biotin-TEG

**Bases Option:** Phosphorothioate bond


**Wobble Bases:** Available at no additional cost

**Format:** dried in tubes

**Free shipping**

Eurogentec all along your NGS process

EUROGENTEC  
**NGS**  
Oligonucleotides



>Note

\*Larger amounts are available on request

# RNAi oligos

## ⇒ WHAT IS RNAi INTERFERENCE?

RNA interference is a mechanism of gene silencing at the mRNA level. This phenomenon is triggered by small interfering (si)RNAs and micro (mi)RNAs.

siRNAs and miRNAs regulate gene expression. They can activate the degradation of the targeted mRNA or prevent its translation. ■

### SPECIFICATIONS

**Length:** From 21 to 27 bases

**Delivered quantity:** 3 nmol • 7 nmol • 12 nmol • 22 nmol • 40 nmol  
• 60 nmol • 80 nmol • 200 nmol\*

**Backbone:** RNA, LNA®, 2'-O-Me RNA, 2'-O-MOE RNA and all linkages

**Modifications:** 5': Phosphate, 6-FAM, Cy<sup>3</sup>, Cy<sup>5</sup>, TET, HEX, ...  
3': DABCYL, TEG-Cholesteryl, TAMRA...

**Purifications:** SePOP Desalting, IEX-RP/HPLC or *in vivo*

**Quality Control:** MALDI-TOF MS


**Format:** Dried

**Packaging:** 2 mL tube

**Documentation:** Technical data sheet

**siRNA Design:** Free and guaranteed

**Shipping:** At room temperature



#### DID YOU KNOW?

Eurogentec can synthesise highly modified and very stable **RNA oligos**.  
Contact us at:  
[unique@eurogentec.com](mailto:unique@eurogentec.com)

Design  
assistance

[info@eurogentec.com](mailto:info@eurogentec.com)

#### >Note

- The antisense strand must either have a free 5'-OH (by default) or 5'-phosphate terminus.
  - Certain modifications can sometimes be useful to increase stability or cellular uptake e.g. Modifying siRNA with cholesterol is used to facilitate tissue / cellular uptake.
  - Various fluorescent dyes can be coupled to the 5'-end of the sense strand oligonucleotide to track transfection efficiency of the corresponding duplex.
- \* 3, 12, 40, 80 nmol for HPLC purification. 7, 22, 60, 200 nmol for SePOP purification. Larger synthesis scales are available on request.

## siRNA

### Custom siRNA Duplexes

Eurogentec has co-developed an exclusive siRNA design platform. PhD-level scientists of our design team use this reliable interface to design custom siRNA for any target of your choice.

Eurogentec **guarantees up to 80% minimum silencing** of your gene of interest with at least one of the 3 duplexes designed and synthesised. ■

### Control siRNA Duplexes

In order to monitor your siRNA experiment conditions, Eurogentec provides siRNA control duplexes and kits including negative and positive controls necessary to validate your experiment.

■ **Negative controls** are siRNA molecules **presenting no homology with any known eukaryotic gene**. siRNA controls are already annealed and shipped lyophilised solution. The sequence is properly validated.

■ **Positive controls** consist of siRNA **directed against a range of endogenous and reporter genes**. They are available in 5nmol final quantities. Each control contains 1 siRNA duplex. All siRNA control duplexes are PAGE purified and 100% MALDI-TOF Mass Spectrometry controlled. The sequences are validated and published. ■

Catalogue control siRNA	
Name	Reference
Control siRNA duplex negative control	SR-CL000-005
Control siRNA duplex LaminB1 (human)	SR-CL001-005
Control siRNA duplex Vimentin (human)	SR-CL002-005
Control siRNA duplex NuMA (human)	SR-CL003-005
Control siRNA duplex Beta-actin (human)	SR-CL004-005
Control siRNA duplex Eg-5 (human)	SR-CL005-005
Control siRNA duplex Cdk-1 (human)	SR-CL006-005
Control siRNA duplex pGL2 luciferase (firefly)	SR-CL010-005
Control siRNA duplex pGL3 luciferase (firefly)	SR-CL011-005
Control siRNA duplex GFP (jellyfish)	SR-CL020-005
Control siRNA duplex <i>in vivo</i> 12 nmol	SR-IVSR-004
Control siRNA duplex <i>in vivo</i> 40 nmol	SR-IVGR-020

**CUSTOM RNA OLIGO  
CAN BE ORDERED ONLINE  
VIA THE CUSTOM OLIGO  
CONFIGURATOR.**

Principle of siRNA-mediated RNA interference. The annealed siRNA enter the cell (1). Once inside, double stranded RNA is recognised by the RISC complex. Sense strand siRNA is displaced and the mRNA anneal to the antisense siRNA fixed to the RISC complex (2). mRNA is digested (3) and the RISC complex containing the siRNA is then recycled to begin a new cycle (4).

## miRNA

miRNA (for microRNA) are natural small non-coding RNAs forming short hairpins. They are implied in gene expression and RNA silencing.

### Clear-MiR™ miRNA Inhibitors


Clear-MiR™ miRNA inhibitors are chemically modified antisense RNA oligonucleotides optimised to specifically **target miRNA molecule in cells**. ■

### Add-MiR™ miRNA Mimics


Add-MiR™ oligonucleotides are custom double-stranded synthetic miRNA **mimicking the action of endogenous miRNAs**. ■

#### DID YOU KNOW

2'-O-Me RNA base and phosphorothioate links bring to the RNA oligo a higher stability and a resistance against nuclease. ■



**>Note**  
Clear-MiR™ miRNA Inhibitors and Add-MiR™ miRNA Mimics are available with different labels and can be linked to cholesterol to increase cellular uptake. On request, peptides can also be covalently linked.



# qPCR Probes

Eurogentec offers a wide range of fluorophores and quenchers in various combinations to fit any method and Real-Time thermocycler. ■

## SPECIFICATIONS

**Length:** From 15 to 50 bases

**Synthesis scale:** 10 nmol • 40 nmol • 200 nmol  
• 1000 nmol • 2.5 µmol • 5 µmol • 10 µmol\*

**Backbone:** DNA, LNA®, 2'-O-Me RNA and phosphodiester linkage

**Modifications:** 5': 6-FAM, HEX, Cy<sup>®</sup>3, TET, Cy5<sup>®</sup>...  
3': TAMRA, DABCYL, BHQ™, DDQ...

**Purifications:** RP-HPLC or Dual HPLC

**Quality Control:** MALDI-TOF MS and analytical HPLC

**Format:** Dried

**Packaging:** 2 mL tube

**Documentation:** Technical data sheet

**Probe Design:** Available on request

**Shipping:** At room temperature

### FOR COMMERCIAL USE

**Access™ Dyes** from Eurogentec are a class of high performance proprietary molecules with IP-friendly terms and conditions. ■

### ALSO AVAILABLE

Eurogentec provides kits & reagents for **qPCR assays** (see Amplification brochure). ■


For **other modified oligonucleotides** please refer to the custom oligonucleotide chapter (p.9). ■

1 PROBE ORDERED  
=  
2 PRIMERS OFFERED<sup>1</sup>

#### >Note

1. Double-Dye or Molecular Beacon Probe. Each primer must be 15-30 DNA bases, unmodified and RP-Cartridge purified. The synthesis scale of the primers must be similar to the probe one (10, 40, 200 or 1000 nmol). Offer Valid only for orders placed on the Eurogentec e-commerce platform EOS. General conditions of sale will be applied.

\* Larger synthesis scales are available on request




## ⇒ DOUBLE-DYE PROBES

**EUROGENTEC OFFERS** a large range of fluorescent dyes and quenchers including proprietary efficient molecules: HiLyte Fluor™ dyes and QXL™ quenchers.

### LNA® Double-Dye probes

LNA® bases have a modification to the ribose backbone that locks the base in the C3'-endo position, which favors RNA A-type helix duplex geometry. Compared to DNA Double-Dye probes, LNA®



Double-Dye probes exhibit higher thermal stabilities, specificity and reproducibility. They show better mismatch discrimination which allows the use of shorter probes.

Furthermore, LNA® offers the possibility to adjust  $T_m$  values of primers and probes in multiplex assays. ■




## ⇒ MOLECULAR BEACONS

**EUROGENTEC IS** a licensed supplier of Molecular Beacons and offers standard, wavelength-shifting and 2' O-Me RNA molecular beacon.

### 2' O-Methyl RNA Molecular Beacons


2' O-Methyl RNA probes perform better than DNA oligonucleotides. They are more nuclease resistant, have a higher affinity, specificity and hybridisation kinetics compared to DNA homologues.



### ALSO AVAILABLE

Plexor™ primers. ■

### LC HYBRIDISATION PROBES




Two LC hybridisation probes labelled with a single fluorescent molecule specifically recognise two adjacent sequences in the target DNA. When the probes are bound to the target sequence, the fluorescent signal is transferred from the donor to the acceptor, which starts to fluoresce. A 3' phosphate group is also added to prevent extension of the reporter probe by Taq DNA polymerase during the PCR cycles. ■

>Note

PNA FISH Probes are also available to detect chromosome aberrations in the centromer. Please see p.19 for more information.




## ⇒ MGB PROBES



**EUROGENTEC PROVIDES** high quality MGB probes perfectly suited for patient management [1]. MGB increases the  $T_m$  of a probe because of its minor groove binding ability. MGB probes are more specific, more efficient and more sensitive than standard double-dye probes.

We provide a complete offer with more than 15 dyes covering all qPCR channels - FAM,

- TET,
- AP5,
- Yakima Yellow<sup>®</sup>,
- Texas Red<sup>®</sup>,
- Cy<sup>®</sup>5,
- ROX,
- DragonFly<sup>™</sup> Orange,
- ATTO
- HEX
- JOE



Our MGB Probes are RP-HPLC purified and can be delivered in 6, 20 and 50 nmol, dried or in solution (TE or H<sub>2</sub>O). For maximal convenience, a 10 nmol dried aliquoted format is also available for the 20 and 50 nmol quantities at no additional cost. The probes are quality controlled by MALDI-TOF MS + HPLC and are available in IVD grade upon request. ■

### Legend

- Fluorophore
- Quencher
- MGB
- Polymerase
- Signal record

The dispensing service, in line with ISO15189, brings to your in-house PCR/qPCR assays or commercial kits, a diagnostic-like grade fitting the highest quality standards. **Any size of routine assays to full kitting solutions** can be produced with a very high reliability, reproducibility and accuracy. This process saves set-up time and reduces reagent wastage, while keeping format flexibility.



**Eurogentec**  
Experience true partnership



## High-throughput Dispensing service

**Contact:**  
[dispensing@eurogentec.com](mailto:dispensing@eurogentec.com)



Food  
GMO  
Clinical diagnostics  
Veterinary pathogen  
detection  
Custom Buffer  
Formulation


### >Note

<sup>1</sup> Restriction of use in the following countries: FR, UK, DE, IT, JP, ES. In these countries MGB probes must only be used for patient management. Use is free of limitation in other countries. End users are covered under Eurogentec's conveyed license for patient management.





**Eurogentec**  
Experience true partnership



Oligos

# DYES & QUENCHERS COMPATIBILITY TABLE

> Highly complex oligonucleotides > Real-Time qPCR


## FLUORESCENT DYES

## QUENCHERS

Research & Commercial use

Research use only

Emission color	Fluorescent Dyes	Max Abs in nm	Max Em in nm
500	Cascade blue	396	410
	Alexa Fluor® 405	401	421
	AMCA-X	353	442
	Alexa Fluor® 350	346	442
	Pacific Blue™ 416	451	459
	Marina blue®	362	459
	Epoch Blue	381	460
	ATTO 390	390	476
	ATTO 425	439	485
	EDANS	336	490
	BODIPY® 493/503	500	504
	ATTO 465	453	509
	CY 2	490	510
	BODIPY® FL	505	513
	Alexa Fluor® 488	495	519
	FAM	495	520
	DY-495	493	521
	ATTO 488	500	520
	Hilyte™ Fluor 488	497	525
	Fluorescein	494	525
	Oregon Green® 488-X	495	525
	Oregon Green® 488	501	526
	Oregon Green® 514	489	526
	ATTO 495	498	526
	Oyster-500	500	527
	Rhodamine Green™	505	527
	Rhodamine Green™-X	503	528
	DY 505	505	530
	Alexa Fluor® 430	433	541
	TET	521	536
	ATTO 520	517	538
	Lucifer Yellow	428	532
	Alexa Fluor® 514	517	542
	CAL Fluor Gold 540	522	544
	Yakima Yellow®	530	549
AP5	527	549	
BODIPY® R66	528	550	
ATTO 532	532	552	
Alexa Fluor® 532	532	553	
BODIPY® 530/550	534	554	
VIC	538	554	
JOE	529	555	
HEX	535	556	
CAL Fluor™ Orange	540	559	
DY-488XL	485	560	
Alexa Fluor® 555	555	565	
Hilyte™ Fluor-555	550	566	
Quasar 570	548	566	
BODIPY® 558/568	558	569	
CY 3	550	570	
Oyster-554	556	570	
BODIPY® 564/570	565	571	
DY 564	551	572	
DY-555	547	572	
DY-556	548	573	
Alexa Fluor® 546	556	573	
BODIPY® TMR	542	574	
NED	546	575	
Dragonfly orange™	554	576	
ATTO 550	554	576	
DY 540	559	578	
Rhodamine Red™-X	560	580	
TAMRA	557	583	
Rhodamine Red™	570	590	



**ORDER**

Oligonucleotides can be ordered via the web oligo configurator or on request at [unique@eurogentec.com](mailto:unique@eurogentec.com). For more information, please contact us at [oligo@eurogentec.com](mailto:oligo@eurogentec.com)

**HIGH PERFORMANCE**

**Access™ Dyes and Quenchers\*** from Eurogentec are a class of high performance proprietary molecules with IP-friendly terms and conditions. ■

**SEE ALSO**

**Compatibility cycler channel vs fluorescent dyes page 27.**


**USAGE**

**In bold blue:** Dyes & Quenchers for commercial & research use in black: Dyes & Quenchers for research use only (RUO)  
 In light grey: Not available at Eurogentec

- Dark blue bands represent the indicative quenching range of recommended quenchers.
- Gray bands represent the indicative quenching range of alternate quenchers.
- Purple bars represent the optimal recommended quenching range.
- ★ Access™ Dyes & Quenchers

All wavelength values are in nm.

Note: Spectral properties of a dye or a quencher may vary when coupled to another molecule



DY-510XL	509
BODIPY™ 576/589	576
<b>ATTO 565</b>	<b>564</b>
CAL Fluor™ Red 590	569
BODIPY™ 581/591	586
FRET 588	595
CY3.5	591
DY-590	590
<b>ROX 575</b>	<b>602</b>
Alexa Fluor™ 568	578
CAL Fluor™ Red 610	590
<b>AP593</b>	<b>593</b>
LC Red 610	610
Texas Red™	597
<b>HiLyte™ Fluor 594</b>	<b>593</b>
BODIPY™ TR	589
Alexa Fluor™ 594	590
<b>ATTO 590</b>	<b>593</b>
ATTO 594	603
Alexa Fluor™ 610	612
DY-610	610
DY-680XL	500
ATTO 610	616
CAL Fluor™ Red 635	618
BODIPY™ 630/650	625
LC Red 640	625
DY-615	621
ATTO 620	620
Alexa Fluor™ 633	632
PULSAR 650	640
ATTO 633	630
DY-630	636
DY-632	637
DY-633	637
DY-631	637
BODIPY™ 650/665	646
DY-601XL	606
DY-520-XL	520
ATTO 647N	646
Alexa Fluor™ 647	650
Oyster-645	645
CY5	650
Quasar™ 670	647
WellRed D4-PA	650
DY-635	647
DY-636	645
DY-647P1	653
DY-650	653
Oyster-656	656
<b>HiLyte™ Fluor 647</b>	<b>650</b>
DY-651	656
ATTO 655	663
Alexa Fluor™ 660	663
CY5.5	675
IRDye 700	680
ATTO 680	681
DY-675	674
DY-676	674
Alexa Fluor™ 680	679
Quasar™ 705	690
WellRED D3-PA	685
DY-681	691
DY-680	690
LC Red 705	680
ATTO 700	700
Alexa Fluor™ 700	702
DY-700	707
DY-701	706
ATTO 725	728
DY-730	732
DY-732	736
DY-731	736
ATTO 740	743
WellRED D2-PA	750
Alexa Fluor™ 750	749
CY7	747
DY-750	747
<b>HiLyte™ Fluor 750</b>	<b>753</b>
DY-751	751
DY-776	771
DY-781	783
DY-782	783
DY-780	770
IRDye 800	787



BECAUSE YOUR EXPERIMENTS REQUIRE ALWAYS MORE CUSTOMISATION, UNIQUE OLIGONUCLEOTIDES BRING YOU THE PERFECT SOLUTION.

# Unique™ oligos

## ⇒ HIGHLY COMPLEX OLIGOS

If you cannot find here the oligonucleotides that best suits your needs please think about unique oligos.

Eurogentec can synthesise highly complex oligonucleotides and can incorporate any of the modifications from recognised chemical suppliers (Glen research, TriLink BioTechnologies...).

Send us the specifications of your Unique™ oligonucleotides (sequence or length, chemistries, modifications, purifications, expected purity, synthesis scale or final amount, format, packaging...) or fill the form online (<https://secure.eurogentec.com/unique-oligonucleotides-quote-request-form.html>) and you will receive the corresponding information in terms of technical feasibility, pricing and turnaround times within 2 working days. ■



### SPECIFICATIONS

**Length:** From 2 to 225 bases

**Synthesis scale:** Customised

**Backbone:** Usual or atypical chemistry

**Modifications:** Common or rare modifications

**Purifications:** SePOP desalting, RP-Cartridge•Gold™, HPLC, PAGE, Dual HPLC, UltraPureGold™

**Quality Control:** Adapted to your needs

**Format:** Adapted to your needs

**Packaging:** Adapted to your needs

**Documentation:** Technical data sheet custom documentation

**Shipping:** As defined by the customer

### ALSO AVAILABLE

#### Custom Gene Synthesis

- From simple gene to highly complex sequence
- Up to 50 kbp
- Gene Optimisation
- 100% Guaranteed sequence
- Fast turnaround time

Contact: [gene@eurogentec.com](mailto:gene@eurogentec.com)

More info in the [SmartGene Brochure](#)

### CONTACT US

[unique@eurogentec.com](mailto:unique@eurogentec.com)

## Oligonucleotides | Aptamers

**Eurogentec and Novaptech** join their expertise to provide high-value aptamer service

- > Good alternative to antibodies
- > Fast and controlled production
- > Any target - various conditions

Ingeniously shaped to fit

# Catalogue

## oligos

Name	Sequence	Bases	T <sub>m</sub> (°C)	Ref.
16S rRNA For	AGA GTT TGA TCC TGG CTC AG	20	55.2	UN-PR001-005
16S rRNA Rev	ACG GCT ACC TTG TTA CGA CTT	21	57.4	UN-PR005-005
3' RACE PCR	GGC CAC GCG TCG ACT AGT AC	20	60.6	UN-PR010-005
Anchored Oligo dT (20)	TTT TTT TTT TTT TTT TV	20	39.2	UN-PR015-005
Anchored Oligo dT (22)	TTT TTT TTT TTT TTT TV N	22	42.8	UN-PR020-005
Bluescript KS	TCG AGG TCG ACG GTA TC	17	53.3	UN-PR025-005
Bluescript SK	CGC TCT AGA ACT AGT GGA TC	20	52.4	UN-PR030-005
cDNA Cloning Primer	GGC CAC GCG TCG ACT AGT ACT TTT TTT TTT TTT TV	38	64.8	UN-PR035-005
EGFP-C	CAT GGT CCT GCT GGA GTT CGT G	22	61.2	UN-PR040-005
EGFP-N	CGT CGC CGT CCA GCT CGA CCA G	22	67.2	UN-PR045-005
G3PDH For	ACC ACA GTC CAT GCC ATC AC	20	58.6	UN-PR050-005
G3PDH Rev	TCC ACC ACC CTG TTG CTG TA	20	59.7	UN-PR055-005
M13 Forward (-20)	GTA AAA CGA CGG CCA GT	17	53.0	UN-PR060-005
M13 Forward (-41)	CGC CAG GGT TTT CCC AGT CAC GAC	24	65.5	UN-PR065-005
M13 Reverse (-27)	CAG GAA ACA GCT ATG AC	17	47.3	UN-PR070-005
M13 Reverse (-48)	AGC GGA TAA CAA TTT CAC ACA GG	23	57.2	UN-PR075-005
Neomycin For	CTT GGG TGG AGA GGC TAT TC	20	55.6	UN-PR080-005
Neomycin Rev	AGG TGA GAT GAC AGG AGA TC	20	54.0	UN-PR085-005
Oligo dT, 15mer	TTT TTT TTT TTT TTT	15	29.7	UN-PR090-005
Oligo dT, 16mer	TTT TTT TTT TTT TTT T	16	32.1	UN-PR095-005
Oligo dT, 18mer	TTT TTT TTT TTT TTT TTT	18	36.0	UN-PR100-005
Oligo dT, 20mer	TTT TTT TTT TTT TTT TTT TT	20	39.1	UN-PR105-005
PCMV Forward	CGC AAA TGG GCG GTA GGC GTG	21	64.8	UN-PR110-005
pET 3'	CTA GTT ATT GCT CAG CGG	18	50.6	UN-PR115-005
pET 5' (T7)	TAA TAC GAC TCA CTA TAG G	19	45.3	UN-PR120-005
pET Upstream	ATG CGT CCG GCG TAG A	16	56.7	UN-PR125-005
pGEX 3'	CCG GGA GCT GCA TGT GTC AGA GG	23	65.2	UN-PR130-005
pGEX 5'	GGG CTG GCA AGC CAC GTT TGG TG	23	67.0	UN-PR135-005
ROSA26 Promoter For	AAA GTC CTG CTT AGT TGT TAT	21	53.2	UN-PR140-005
ROSA26 Promoter Rev	GGA GCG GGA GAA ATG GAT ATG	21	56.3	UN-PR145-005
SP6 Promoter	TAC GAT TTA GGT GAC ACT ATA G	22	50.0	UN-PR150-005
SP6 Upstream	ATT TAG GTG ACA CTA TAG	18	42.8	UN-PR155-005
T3 Promoter	AAT TAA CCC TCA CTA AAG GG	20	50.4	UN-PR160-005
T7 Promoter	TAA TAC GAC TCA CTA TAG GG	20	48.3	UN-PR165-005
T7 Terminator	GCT AGT TAT TGC TCA GCG G	19	54.1	UN-PR170-005

## ⇒ UNIVERSAL PRIMERS

Universal primers are complementary to nucleotide sequences that occur very commonly in specific sets of DNA molecules and cloning vectors. Thus, they are able to **bind to a wide variety of DNA templates.** ■

## SPECIFICATIONS

**Quantity:** 1 OD/5 nmol

**Backbone:** DNA

**Modifications:** None

**Purifications:** RP-HPLC

**Quality Control:** MALDI-TOF MS + CGE

**Format:** Dried

**Packaging:** 2 mL tube

**Documentation:** Technical data sheet

**Shipping:** At room temperature

## ⇒ PNA FISH

In principle, fluorescence *in situ* hybridisation (FISH) should be able to provide information on the telomere length of individual chromosomes. The efficiency of conventional labelled oligos is not sufficient to be extended beyond qualitative studies of TTAGGG repetitions in chromosomes of various species.

PNA chemical structure brings a **higher sequence specificity**, an improved **stability**, better **reproducibility**, and lower background noise. Due to the higher T<sub>m</sub> of PNA/DNA duplexes, short (18-mer) telomere PNA (CCCTAA)<sub>3</sub> are now widely used. ■

### SPECIFICATIONS

**Length:** 18 bases

**Quantity:** 5 nmol

**Backbone:** PNA

**Modifications:** FAM • Cy3<sup>®</sup> • Cy5<sup>®</sup> • FITC • TMR

## ⇒ CALIBRATION OLIGOS

**Dye-labelled calibration oligos** are a set of 5' fluorescent dT10 oligonucleotides recommended to calibrate some real-time qPCR thermocyclers. Calibration is a preliminary step indicated to adjust fluorescent signal analysis. qPCR Dye Calibration Oligos' enables the thermocycler to recognise the spectra of each single dye and to control signal overlap that may occur in multiplexed assays particularly. ■

### SPECIFICATIONS

**Length:** 10 bases

**Quantity:** 5 nmol

**Backbone:** DNA

**Modifications:** AP5, Yakima Yellow<sup>®</sup>, HEX, Dragonfly Orange<sup>™</sup>, TET, JOE, HiLyte<sup>™</sup> Fluor 647, ROX

**Purification:** RP-HPLC

**Quality control:** MALDI-TOF MS

**Format:** Dried

**Packaging:** 2 mL tube

**Documentation:** Technical data sheet

**Shipping:** At room temperature



### ORDERING INFORMATION

Name	Quantity	#Cat	Name	Quantity	#Cat
<b>C-Rich Telomere Probes</b>			<b>Centromere Probes</b>		
TelC-FAM	5 nmole	PN-TC001-005	Cent-Cy3	5 nmole	PN-CN050-005
TelC-Cy3	5 nmole	PN-TC050-005	Cent-FAM	5 nmole	PN-CN001-005
TelC-Cy5	5 nmole	PN-TC055-005	Cent-Cy5	5 nmole	PN-CN055-005
TelC-Alexa488	5 nmole	PN-TC060-005	Cent-Alexa488	5 nmole	PN-CN060-005
TelC-FITC	5 nmole	PN-TC011-005	Cent-FITC	5 nmole	PN-CN011-005
TelC-TAMRA	5 nmole	PN-TC030-005	Cent-TAMRA	5 nmole	PN-CN030-005
TelC-Alexa647	5 nmole	PN-TC020-005	Cent-Alexa647	5 nmole	PN-CN020-005
TelC-Biotin	5 nmole	PN-TC040-005	Cent-Biotin	5 nmole	PN-CN040-005
<b>G-Rich Telomere Probes</b>			<b>Centromere Protein B Probes</b>		
TelG-FAM	5 nmole	PN-TG001-005	CENPB-FAM	5 nmole	PN-CP030-005
TelG-Cy3	5 nmole	PN-TG050-005	CENPB-Cy3	5 nmole	PN-CP050-005
TelG-Cy5	5 nmole	PN-TG055-005	CENPB-Cy5	5 nmole	PN-CP055-005
TelG-Alexa488	5 nmole	PN-TG060-005	CENPB-Alexa488	5 nmole	PN-CP060-005
TelG-FITC	5 nmole	PN-TG011-005	CENPB-FITC	5 nmole	PN-CP011-005
TelG-TAMRA	5 nmole	PN-TG030-005	CENPB-TAMRA	5 nmole	PN-CP001-005
TelG-Alexa647	5 nmole	PN-TG020-005	CENPB-Alexa647	5 nmole	PN-CP020-005
TelG-Biotin	5 nmole	PN-TG040-005	CENPB-Biotin	5 nmole	PN-CP040-005



### ORDERING INFORMATION

Name	Sequence	Modification 5'	Bases	Abs/Em (nm)	Reference
AP5-T10	TTTTTTTTTT	AP5	10	527/549	UN-CT001-005
YY-T10	TTTTTTTTTT	YY	10	530/549	UN-CT005-005
HEX-T10	TTTTTTTTTT	HEX	10	535/556	UN-CT010-005
DFO-T10	TTTTTTTTTT	DFO	10	554/576	UN-CT015-005
TET-T10	TTTTTTTTTT	TET	10	521/536	UN-CT020-005
JOE-T10	TTTTTTTTTT	JOE	10	529/555	UN-CT025-005
HL647-T10	TTTTTTTTTT	HL647	10	650/675	UN-CT030-005
ROX-T10	TTTTTTTTTT	ROX	10	575/602	UN-CT035-005

# Glen Synthesis reagents

Eurogentec is an **authorised distributor for Glen research** products in most of European countries. We distribute a large range of reagents for the DNA/RNA oligonucleotide synthesis including phosphoramidites, solid supports for oligonucleotide purifications, labelling dyes and modifications. ■

- CE Phosphoramidites, synthesis columns and solvents/reagents
- Other monomers
- Minor bases
- Modification and labelling
- RNA Synthesis
- And many more

 **GLEN**  
RESEARCH  
part of Maravai LifeSciences



# Additional Services

## Additional QC

**MALDI-TOF Mass Spectrometry:** This method provides the most precise information about the length, deprotection-product and the presence of labels for modified oligonucleotides over a broad range of lengths (up to 60 bases).

**RP-UHPLC:** This is a very efficient technique giving quantitative information about the purity level of oligonucleotides from 15 to 40 bases long.

**IEX-UHPLC:** This technique is particularly adapted to quantify the purity level of oligonucleotides from 15 to 40 bases long.

**Capillary Gel Electrophoresis (CGE):** This method is adapted to assess very precisely the purity of oligonucleotides longer than 40 bases (on request).

**Fluorescence analysis:** This non-destructive physical technique provides qualitative information about your fluorescent oligonucleotides. ■

## Format

**Dried:** All the synthesised oligonucleotides are dried by default (except SePOP unmodified oligonucleotides from 15 to 39 bases).

**In solution:** You may select the nature of the reconstitution buffer (H<sub>2</sub>O or TE), the volume of the reconstitution buffer (from 50 to 1000 µl) or/and the final oligonucleotides concentration (from 5 to 250 µM).

**Annealed:** siRNA or cloning linkers are annealed by default.

**Mixed:** Similar amounts of forward and reverse oligonucleotides can be mixed in a single tube. ■

## Packaging

**2 mL tube:** By default, each oligonucleotide is provided in individual 2 mL tube. Higher volume can be delivered on request (15 mL, 50 mL)

**96-well plates:** Cluster tubes, well plate and deep well plate are available.

**384-well plates:** Specially suitable for high throughput experiments requiring more than 96 oligonucleotides.

**Aliquoting:** All the oligonucleotides in solution can be split in small aliquots of desired volume (from 50 to 1000 µl). ■

## Shipping

Your oligonucleotides can be express shipped in 24 hours upon request (see page 24 for more details). ■

## Design

We continuously update our software and design rules to reflect the latest scientific developments as well as integrate customer requirements. This service includes primers, Double-Dye Oligonucleotides, Molecular Beacons, siRNA design, miRNA inhibitors... ■

## Increase the quality of your oligo for demanding applications

### Hospital and commercial kits

#### Track™ oligonucleotides

Track™ oligonucleotides offer a higher traceability (and other quality assets) in the production process than life science research oligonucleotides.

#### cGMP oligonucleotides

cGMP oligonucleotides ensure exceptional product quality by manufacturing in classified cleanrooms and use of an ISO 13485-certified and GMP compliant QMS.

### Therapeutics fields

#### Pre-clinical oligonucleotides

Large scale pre-clinical oligonucleotides are manufactured in cleanrooms and delivered with appropriate documentation. Additional QC tests such as endotoxin level are offered. For more information download our IVD brochure.

[www.eurogentec.com/invitrodiagnostics.html](http://www.eurogentec.com/invitrodiagnostics.html)

OLIGO GRADE					
Lab Services Hospital & Clinical Labs	Discovery	Routine Assays		Contamination Sensitive Assays	
Diagnostic Companies	Discovery	Feasibility	Prototyping	Validation	Commercialisation
Oligonucleotide Grades	Research	Track	Pre-Diagnostic	Diagnostic	
<b>Process</b>					
Dedicated Account Contact Person	Option	✓	✓	✓	✓
Customised Fill & Finish	Option	Option	✓	✓	✓
<b>Quality Management</b>					
ISO 9001 Certification	✓	✓	✓	✓	✓
ISO 13485 Certification	-	-	✓	✓	✓
Qualification/Validation [Equipment & Method]	-	Partial	✓	✓	✓
<b>Control</b>					
Quantification	Single	Dual	Triple	Triple	Triple
Stringent QC Tests (validated)	-	✓	✓	✓	✓
Traceability	Partial	Documented	Documented	Documented	Full documented
Batch Record [Archived for 5 years]	-	-	Partial	Full	Full
Classified Cleanroom	-	-	✓	✓	✓
Certificate of Analysis [CoA]	-	✓	✓	✓	✓



## ANNEXES

## SYNTHESIS SCALE VS GUARANTEED YIELD

## GOOD TO KNOW

All allowed purifications are represented in this table. To select the recommended purification according to your applications and modifications, please refer to p.8.

The **synthesis scale** refers to the **amount of raw material** used to start the synthesis of oligonucleotides. The **yield** corresponds to the amount of **final product**

recovered at the end of the synthesis and purification processes. The length, the sequence, the type/number of modifications and the purification, strongly

influence the reaction yield. Based on that, Eurogentec defined a minimum guaranteed yield in nmoles for all product categories (see table below). The minimum guaranteed yields

represent only a reference because the delivered quantities may vary.

		Synthesis scale (nmol)														Minimum Guaranteed Yield (nmol)																					
		10	40	200	1000	2500	5000	10000	20000	Purification																											
Range	Product	Length	SePOP	RP-Cartridge-Gold™ HPLC (RP or IEX)	RP-Cartridge-Gold™ HPLC (RP or IEX)	PAGE <sup>(3)</sup>	Dual HPLC	SePOP	RP-Cartridge-Gold™ HPLC (RP or IEX), in vivo	PAGE <sup>(3)</sup>	Dual HPLC	UltraPureGold™	SePOP	RP-Cartridge-Gold™ HPLC (RP or IEX), in vivo	PAGE <sup>(3)</sup>		Dual HPLC	UltraPureGold™	SePOP	RP-Cartridge-Gold™ HPLC (RP or IEX)	Dual HPLC	PAGE <sup>(3)</sup>	SePOP	RP-Cartridge-Gold™ HPLC (RP or IEX)	Dual HPLC	PAGE <sup>(3)</sup>	SePOP	RP-Cartridge-Gold™ HPLC (RP or IEX)	Dual HPLC	PAGE <sup>(3)</sup>	SePOP	RP-Cartridge-Gold™ HPLC (RP or IEX)	Dual HPLC	PAGE <sup>(3)</sup>			
Custom Oligonucleotides	Non-Modified (DNA only)	5-9	-	-	-	-	-	60	50	30	20	15	-	180	100	80	40	40	-	450	200	100	100	1000	400	200	200	1800	800	400	400	-	-	-	-		
		10-19	5	4	-	20	16	10	4	3	70	60	45	30	23	15	200	140	100	70	50	30	500	250	125	125	1000	500	250	250	2000	1000	500	500	4200	2100	1050
		20-39	5	4	-	20	16	10	4	2	60	50	30	20	15	10	190	120	90	40	45	20	475	225	115	115	1000	500	250	250	2000	1000	500	500	4200	2100	1050
		40-59	3	2	-	10	8	5	2	1	30	25	15	12	7	6	115	60	45	20	20	12	285	110	55	55	600	230	115	115	1200	460	230	230	2500	1000	500
		60-79	2	2	-	8	6	-	2	-	20	18	-	8	-	4	75	40	-	14	-	8	185	-	-	40	350	-	-	90	750	-	-	180	1500	-	-
		80-99	-	-	-	-	-	-	1	-	-	-	-	3	-	2	-	-	-	5	-	3	-	-	-	30	-	-	-	40	-	-	80	-	-	-	
100-139	-	-	-	-	-	-	-	-	-	-	-	2	-	1	-	-	-	3	-	2	-	-	-	10	-	-	-	20	-	-	40	-	-	-			
Custom Oligonucleotides	Modified <sup>(1)</sup> (including DNA, RNA, 2' O-Me RNA, LNA and phosphorothioate linkages)	5-9	-	-	-	3	-	-	-	12	-	6	-	-	25	-	12	-	-	60	30	30	-	125	60	60	-	250	125	125	-	-	-	-			
		10-19	-	-	-	12	6	5	4	1	35	20	17	15	8	-	70	40	35	30	15	-	175	90	45	45	500	190	95	95	1000	380	190	190	2000	-	-
		20-59	-	-	-	8	5	4	3	1	20	15	12	10	6	-	45	35	25	20	12	-	100	65	30	30	300	135	65	65	600	275	130	130	1200	600	275
		60-139	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	15	-	-	30	-	-	60	-	-	-	-	-	-
Real-Time qPCR Probes	Double-Dye probes <sup>(2)</sup>	8-38	-	-	<2 <sup>(4)</sup>	-	4	-	-	-	12	-	-	-	-	25	-	-	-	-	65	-	-	-	135	-	-	-	275	-	-	-	600	-	-		
	Molecular Beacons	32-50	-	-	1	-	-	-	-	4	-	-	-	-	-	12	-	-	-	-	30	-	-	-	65	-	-	-	130	-	-	-	275	-	-		
	MGB Taqman Probes	8-30	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	30	-	-	-	65	-	-	-	130	-	-	-	275	-	-		
		Delivered quantity: 6, 20 or 50 nmol																																			
RNAi Oligonucleotides	siRNA Duplexes Non-Modified <sup>(5)</sup>	21-27	7	-	3	22	-	12	-	-	60	-	40	-	-	200	-	80	-	-	-	On Request															
	siRNA Duplexes Modified <sup>(1)</sup>	21-27	7	-	3	22	-	12	-	-	60	-	40	-	-	200	-	80	-	-	-	On Request															
NGS Oligonucleotides	RP-Cartridge purified	20-85	Minimum delivered quantity: 10 nmol																																		
	RP-HPLC purified																																				
Universal Primers	-	15-38	-	-	-	5	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-			
Unique Oligonucleotides	-	2-225	On request - please contact us at <a href="mailto:unique@eurogentec.com">unique@eurogentec.com</a>																																		

## Post-synthesis modifications may yield 50% less than the above stated values.

Table: (1) Between 5 and 59 bases length single-modified Oligonucleotides. Eurogentec does not provide minimum guaranteed yield for modified oligonucleotides greater than 59 bases. Post-synthesis modifications are not compatible with SePOP and RP-Cartridge-Gold™ purification. A lower yield may result from poly-modifications and/or strong secondary structures.

(2) Double-Dye probes only result from the combination of a 5' fluorescent dye and a 3' quencher.

(3) Except for oligonucleotides with GC-rich regions.

(4) Only available for Double-Dye FAM-TAMRA 10 nmol and FAM-BHQ1 10 nmol.


(5) Non-modified siRNA's only include 3' dTdT overhang.

## List of the post-synthesis modifications

> 5' Alexa Fluor\* (350, 430, 488, 500, 514, 532, 546, 555, 568, 594, 610, 633, 647, 660, 680, 700 and 750)  
> 5' ATTO (390, 425, 465, 488, 495, 520, 532, 550, 565, 590, 594, 610, 620, 633, 635, 647N, 655, 680, 700, 725 and 740)  
> 5' BODIPY\* (530/550, FL and TR)  
> 3', 5' and dT Cascade Blue\*  
> 3' and dT Cy\* (3, 3.5, 5 and 5.5)  
> 3', 5', dR and dT Digoxigenin  
> 5' Dragonfly Orange\*  
> 5' DY\* (681, 781 and 782)  
> dR 6-FAM

> dR and dT HEX  
> 5' HiLyte™ Fluor (405, 488, 555, 594, 647, 680 and 750)  
> 3', dR and dT JOE  
> 5' Marina Blue\*  
> 5' Oregon Green\* (488 and 488 X)  
> 5' Pacific Blue™  
> 3' QXL\*  
> 3', 5', dR and dT Rhodamine 6G  
> 3', 5', dR and dT ROX  
> 5' TAMRA  
> dR and dT TET  
> 3', 5', dR and dT Texas Red\*

For more information, please contact our Oligo Centre at: [oligo@eurogentec.com](mailto:oligo@eurogentec.com) or visit our website: [www.eurogentec.com](http://www.eurogentec.com)



## FOR PCR

## DETERMINE THE RIGHT SYNTHESIS SCALE

Final Oligo Concentration	50 nM	150 nM	300 nM	600 nM	900 nM	
	Average number of Reactions (total volume 100µL)					Minimum quantity to order*
	100	30	15	8	6	0.5 nmol
	1000	300	165	80	55	5 nmol
	5000	1650	830	415	275	25 nmol
	10 000	3300	1660	830	555	50 nmol
	100 000	33 300	16 660	8330	5555	500 nmol

\*Please select in the minimum guaranteed yield table the synthesis scale corresponding to the desired minimum quantity.

## CHECK YOUR SHIPPING METHODS

The delivery time depends on the specifications of your oligonucleotides (see table below).

### Eco-Logik Delivery

■ By local Mail to reduce the global ecological impact. Receive your oligonucleotides in your mailbox. Available for Belgium, France and Monaco

### Express Delivery

■ All oligonucleotides  
■ By Express courier to receive your oligonucleotides as fast as possible (24 to 48 hours) in your hands.

■ Same day shipping option  
- For orders received before 10.00 AM (Central European Time)  
- For custom oligonucleotides (max 24), 10/40 nmol scale,

5-30 DNA bases, unmodified, SePOP desalted or RP-Cartridge purified.

## Delivery times (in working day)

Range	Product	Length	Purification					
			SePOP	RP-Cartridge-Gold™	HPLC (RP or IEX)	PAGE	Dual HPLC	UltraPureGold™
Custom Oligonucleotides	Non-Modified (DNA Only)	5-9	-	4-5	5	6	7	7
		10-39	2-3	4-5	5	6	7	7
		40-59	5	6	7	8	-	9
		60-79	-	6	-	8	-	9
		80-139	-	-	-	10	-	11
	Modified (including DNA, RNA, 2' O-Me RNA & LNA <sup>(1)</sup> )	10-39	5	7	7	8	9	9
		40-59	7-8	9-10	9-10	10-11	11-12	11-12
Real-Time qPCR Probes	Double-Dye Probes	8-38	-	-	7	-	-	-
	Molecular Beacons	32-50	-	-	12-15	-	-	-
	MGB Taqman Probes	8-30	-	-	5-7	-	-	-
RNAi Oligonucleotides	siRNA Duplexes	21-27	5-7	-	9-10	10-11	-	-
NGS Oligonucleotides	-	20-85	-	4-6	5-7	-	-	-
Universal Primers	-	15-38	-	-	2-3	-	-	-
Unique Oligonucleotides	-	2-225	On Request					

For large order or Unique Oligonucleotides, please feel free to contact us at [oligo@eurogentec.com](mailto:oligo@eurogentec.com) to receive more details in terms of delivery schedules. 5'AP, BSA, HRP or SBP Conjugation : 3-5 WD Extra. Additional Purification or Services: 2 WD Extra ; Fax Ordering: 1 WD Extra

## RECEIVE YOUR DOCUMENTATION

Each oligonucleotide is provided with a technical data sheet. Other documentations could be added depending on the oligonucleotide type. All the documents are sent as pdf files to your shipping email address.

	TDS	MS <sup>(1)</sup>	UHPLC	CGE <sup>(3)</sup>
Custom Oligonucleotides	Unmodified	✓		
	Modified	✓	✓ <sup>(4)</sup>	
	UltraPureGold™	✓	✓	✓
Real-Time qPCR Probes	✓	✓	✓	
RNAi Oligonucleotides	✓	✓ <sup>(4)</sup>		
Universal Primers	✓	✓		
Unique Oligonucleotides	✓	✓	✓ <sup>(5)</sup>	✓ <sup>(5)</sup>
NGS oligos	✓	✓		
Calibration oligos	✓	✓		

TDS: Technical Data Sheet; MS: Mass Spectrometry; HPLC: High Performance Liquid Chromatography; Ultra Performance Liquid Chromatography; CGE: Capillary Gel Electrophoresis.

(1) Always provided up to 60 bases long Oligonucleotides.

(2) If applicable.

(3) Can be substituted by another analytical QC

(4) Except for SePOP desalted oligonucleotides.

(5) Optional.

For technical reasons this general rule may be adapted to provide you with the most suitable and useful documentation.

# HOW TO STORE YOUR OLIGO

## Handling information

Products	Format	Storage	Stability**
Custom Oligonucleotides	Dried	RT	18 months
Real-Time qPCR Probes	TE Buffer (pH 8) or dH <sub>2</sub> O	-20 °C	24 months
RNAi Oligonucleotides	Dried	RT	18 months
Catalogue Primers	TE Buffer (pH 8)* or dH <sub>2</sub> O	-20 °C	24 months
PNA FISH Probes / Custom PNA	Dried	RT	18 months

\* Except for Cy<sup>5</sup> dye labelled oligonucleotides (pH7)

\*\* Please protect from light and avoid freeze/thaw cycles.

Please note that depending on sequence and modifications, the stability of the oligos may vary substantially versus the values given above, which should therefore be considered as indicative.

## RECONSTITUTE YOUR OLIGO

1. **Spin** the tube briefly to collect the pellet in the bottom of the tube.
2. **Add** an appropriate volume of recommended buffer.
3. **Allow** the tube to stand a few minutes.
4. **Vortex** the tube for 15 seconds and spin briefly.
5. **Refer to** the dedicated technical data sheet for more information.



## QUANTIFY YOUR OLIGO

To quantify your oligonucleotides, make an aliquot of the resuspended oligonucleotides to a final volume of 1 mL of dH<sub>2</sub>O and vortex for a few seconds. Measure the absorbance of this dilution at 260 nm ( $A_{260}$ ). Use the formula below to calculate the concentration of oligonucleotides in your stock solution. This formula is valid for an absorption of  $A_{260} < 1.2$ .

Concentration in  $\mu\text{g}/\text{mL} = A_{260} \times \text{dilution factor} \times \text{Weight per OD of stock solution (in } \mu\text{g} / \text{OD)}$ .

1 OD<sub>260</sub> (Optical Density) unit is defined as the amount of oligonucleotide which, when dissolved in a volume of 1.0 mL, results in an absorbance of 1.0 when measured at 260 nm in a 1 cm path-length quartz cuvette. 1 OD<sub>260</sub> unit corresponds to approximately 33  $\mu\text{g}$  of single strand DNA. These relationships, however, can be inaccurate for short fragments of DNA, such as oligonucleotides. Base composition and even linear sequence will affect optical absorbance. Hence the precise value of the OD to mass relationship is unique for each oligonucleotide.

### MEASURE

1.0 OD<sub>260</sub> of CCCCCCCCCC (10 bases) equals 39  $\mu\text{g}$   
 whereas 1.0 OD<sub>260</sub> of AAAAAAAAAA (10 bases) equals only 20  $\mu\text{g}$ .

We carefully measure the OD value for your custom oligonucleotide by measuring the absorption at 260 nm using UV spectrophotometer. This information is provided on the oligonucleotide technical data sheet as the number of OD<sub>260</sub> units. The amount of oligonucleotide expressed in nanomoles and micrograms is derived from the OD measurement.

### CALCULATE

Calculate the number of nanomoles present given an OD reading and extinction coefficient:

$$\text{Nanomoles} = (\text{OD}_{260} / \epsilon_{260}) \times 10^6$$

Example:

1 OD<sub>260</sub> unit of primer M13 Forward,  
 5'-GTA AAA CGA CGG CCA GTG-3'  
 Molar extinction coefficient ( $\epsilon_{260}$ ) = 182,800 L / (mole x cm)  
 Nanomoles =  $(1.0 / 182,800) \times 10^6 = 5.47 \text{ nmoles}$

### CONVERT

Convert the amount in nanomoles to micrograms:

$$\text{Micrograms} = \text{Molecular Weight} \times \text{Nanomoles} \times 10^{-3}$$

Example:

1 OD<sub>260</sub> unit of primer M13 Forward,  
 5'-GTA AAA CGA CGG CCA GTG-3'  
 Molecular Weight = 5558.7  
 Micrograms =  $5558.7 \times 5.47 \times 10^{-3} = 30.4 \mu\text{g}$

## CALCULATE THE DATA

### THE MOLAR EXTINCTION COEF.

$$\epsilon_{260} = 2 \times \left( \sum_1^{n-1} \epsilon_{\text{Nearest Neighbour}} \right) - \sum_2^{n-1} \epsilon_{\text{Individual}} + \sum_1^n \epsilon_{\text{Modification}}$$

where  $\epsilon_{\text{Nearest Neighbour}}$  is the nearest neighbour constant for a pair of bases,  $\epsilon_{\text{Individual}}$  is the constant for an individual base, and  $n$  is the length of the oligonucleotide.

### THE MOLECULAR WEIGHT

$$\text{Anhydrous MW (g/mol)} = \sum_{\text{Individual Base}} \text{MW} + \sum_{\text{Individual Mods}} \text{MW} - 63.98 + 2.016$$

#### For DNA bases:

MW dA = 313.21; MW dC = 289.18; MW dG = 329.21; MW dT = 304.20; MW dU = 290.17; MW dI = 314.19

#### For RNA bases:

MW DNA counterpart + 16.  
When determining the weight of Uracil (rU) start with dU and not dT


#### For LNA bases:

MW DNA counterpart + 16 (+42 for dC)

#### For 2' O-Methyl bases:

MW DNA counterpart + 30.03.  
When determining the weight of mU start with dU and not dT

#### For phosphorothioated bases: MW DNA counterpart + 16.06



## WRITE YOUR SEQUENCE IN YOUR ORDER

### IUB CODE

ACGTA = DNA  
(ACGUA) = RNA  
[ACGUA] = 2' O-Me RNA  
<ACGTA> = 2' O-MOE  
{ACGTA} = LNA®  
A\*C\*(G\*U)\*A = Phosphorothioate links

Mixed bases (also known as degenerate or wobble bases) follow the IUB codes:

D=A/G/T  
M=A/C  
H=A/C/T  
I = Inosine = Universal base  
W=A/T  
R=A/G  
Y=C/T  
V=A/C/G  
S=C/G  
K=G/T  
N=A/G/C/T  
B=C/G/T

Oligonucleotide synthesised with mixed bases gives a final product that is a heterogeneous population of distinct species. MW, T<sub>m</sub> and extinction coefficient may be strongly affected by mixed base addition. Rather than reporting the various values for each component, a single value is given.



>Note

ε modification is not known for all modifications.

# COMPATIBILITY CYCLER CHANNELS VS FLUORESCENT DYES

Thermocycler	Blue channel	Channel 1	Channel 2	Channel 3	Channel 4	Channel 5	Channel 6	Channel 7
GeneAmp SDS 5700		FAM	AP5/YY/DE/VC/TET*	DFO*/TAMRA/NEO				
ABI Prism SDS 7000		FAM	AP5/YY/DE/VC/TET*	DFO*/TAMRA/NEO	ROX	-		
ABI Prism SDS 7000		FAM	AP5/YY/DE/VC/TET*	DFO*/TAMRA/NEO	ROX	-		
ABI Prism SDS 7900HT		FAM	AP5/YY/DE/VC/TET*	DFO*/TAMRA/NEO	ROX	-		
ABI Prism SDS 7300		FAM	AP5/YY/DE/VC/TET*	DFO*/TAMRA/NEO	ROX	-		
ABI Prism SDS 7500		FAM	AP5/YY/HEX*/JOE*/VIC/TET*	DFO*/TAMRA/Cy 3/NEO	ROX/TR	Cy 5/HL647	-	
QuantStudio 3, 5, 6, 7		FAM	AP5/YY/HEX*/JOE*/VIC/TET*	DFO*/TAMRA/Cy 3/NEO	ROX/TR	Cy 5/HL647	-	
QuantStudio 7		FAM	AP5/YY/HEX*/JOE*/VIC/TET*	DFO*/TAMRA/Cy 3/NEO	ROX/TR	Cy 5/HL647	Cy 5.5/ATTO 700	
Vii7		FAM	AP5/YY/HEX*/JOE*/VIC/TET*	DFO*/TAMRA/NEO	ROX/TR*	LIZ/ATTO 633	ATTO 680/Alexa Fluor 680	
StepOne		FAM	AP5/YY/HEX*/JOE*/VIC/TET*	ROX	-	-	-	
StepOnePlus		FAM	AP5/YY/HEX*/JOE*/VIC/TET*	DFO*/TAMRA/Cy 3/NEO	ROX	-	-	
Cycler IQ		FAM	AP5/YY/HEX*/CY 3/TET	DFO*/TAMRA/Cy 3/NEO	ROX/TR	Cy 5	-	
My IQ		FAM	-	-	-	-	-	
i5		FAM	AP5/YY/HEX/JOE/TET	DFO*/TAMRA/Cy 3/NEO	ROX/TR	Cy 5	-	
CFX96		FAM	AP5/YY/HEX/JOE/TET	ROX / TR	Cy 5	Cy 5.5	-	
MiniOpticon		FAM	AP5/YY/HEX*/TET*	-	-	-	-	
DNA Engine Opticon 1		FAM	-	-	-	-	-	
DNA Engine Opticon 2		FAM	AP5/YY/HEX*/TAMRA*/VIC/TET	-	-	-	-	
Chromo 4		FAM	AP5/YY/HEX*/TAMRA*/Cy 3/JOE*/VIC/TET	ROX/TR	Cy 5	-	-	
Mx3000P (choice of 4 filters)		FAM	TET	YY*/HEX/JOE/VC	Cy 3	DFO*/TAMRA/NEO	TR/ROX	Cy 5
Mx3000P* (choice of 5 filters)	Epoch Blue / Alexa Fluor 350	FAM	TET	YY*/HEX/JOE/VC	Cy 3	DFO*/TAMRA/NEO	TR/ROX	Cy 5
Mx4000 (choice of 4 filters)		FAM	TET	YY*/HEX/JOE/VC	Cy 3	DFO*/TAMRA/NEO	TR/ROX	Cy 5
Axiom (choice of up to 6 filters)	ATTO 425	FAM	HEX	Cy 3	ROX	Cy 5	-	
Mastecycler epireplex2		FAM	AP5/YY/HEX/JOE/VIC/TET	-	-	-	-	
Mastecycler epireplex4		FAM	AP5/YY/HEX/TET/JOE/VIC	DFO*/TAMRA	ROX	-	-	
LightCycler 15		FAM	AP5/YY/HEX/JOE/VIC/TET	Cy 5	-	-	-	
LightCycler 2		FAM	AP5/HEX/YY*/JOE*/VIC	TR/LC Red 610	ATTO 620/LC Red 640	Cy 5/LC Red 670	LC Red 705/ATTO 680	
LightCycler 480 I	ATTO 425	FAM	AP5/YY/HEX/JOE/VC	ROX/TR/LC Red 610	ATTO 620/LC Red 640	Cy 5	-	
LightCycler 480 II	ATTO 425	FAM	AP5/YY/HEX/JOE/VC	ROX / TR/LC Red 610	ATTO 620/LC Red 640	Cy 5 / Cy 5.5	-	
LightCycler 96		FAM	AP5/YY/HEX/VC	TR/LC Red 610	Cy 5	-	-	
SmartCycler1		FAM	AP5/YY/Cy 3/JOE/VIC/TET	TR	Cy 5	-	-	
SmartCycler2		FAM	AP5/Cy 3/YY*/JOE*/TET	TR/ROX	Cy 5	-	-	
Rotor-Gene 2000 / 3000		FAM	AP5/YY*/JOE/VIC/TET	TAMRA/ROX/Cy 3.5/TR	Cy 5	-	-	
Rotor-Gene 6000	Epoch Blue / Alexa Fluor 350	FAM	AP5/YY/HEX/JOE/VIC/TET	TAMRA/ROX/Cy 3.5/TR	Cy 5	ATTO 680	-	
PikReal		FAM	AP5/HEX/YY	ROX/TR	Cy 5	-	-	
Otover		FAM	AP5/YY/HEX/JOE/VC	DFO*/TAMRA/NEO	ROX/TR/Cy 3.5*	Cy 5	Cy 5.5	

\* perform a dye calibration for optimal results  
 For complementary information, please refer to instrument manufacturer technical guide or contact us at [scientific.support@eurogentec.com](mailto:scientific.support@eurogentec.com)  
 YY = Yakima Yellow  
 DFO = Dragonfly Orange™  
 TR = Texas Red  
 In blue = Recommended by Eurogentec  
 In grey = Not available at Eurogentec



# HOW TO ORDER

ON LINE

WWW.EUROGENTEC.COM

## SINGLE ORDER

1. Connect to [www.eurogentec.com](http://www.eurogentec.com)
2. Click on the oligonucleotide tab of the order centre screen
3. Select the oligonucleotide type (Custom, Probes, RNAi...)
4. Fill the configurator with your oligonucleotide specifications
5. Add your oligonucleotide into your cart and finalise your order

**Custom Oligonucleotides**

Name (automatically generated if this field is empty)

Sequence (3 bases)

5'  3'

DNA = ACGT; RNA = (ACGU); 2'-O-Me RNA = (ACGU); LNA = (ACGT); Phosphorothioates = A\*(G\*(U\*))

If applicable, please insert 1 in the sequence to indicate the location of the internal modification 1 and 2 for the internal modification 2.

Approximate scale

40 nmol (10-95 bases)  Please select one

5' modification  3' modification

Internal modification 1  Internal modification 2

(Indicates with 1: AAN(AAU)) (Indicates with 2: AA(A)AAU)

Additional DC  Frame

None  Dried

**Important notice: 1. Name & 2. Price Primitives**

For each Probe (Double-Dye or Molecular Beacon) in your cart, you can receive 2 free primers. Each primer must be 15-33 DNA bases, unmodified and RP-Cartridge purified. If you do not have a compatible probe in your cart, please configure it via the Double-Dye Probe or the Molecular Beacon configurator.

Please select the Name of the Probe included in these criteria:

## MULTIPLE/ BATCH ORDERS

1. Connect to [www.eurogentec.com](http://www.eurogentec.com)
2. Click on the oligonucleotide tab of the order centre screen
3. Select the Multiple/Batch Order
4. Download the Excel File and fill it in
5. Upload the completed file on the Eurogentec website: [www.eurogentec.com](http://www.eurogentec.com)

**Multiple / Batch Order Oligonucleotides**

1. Download this Excel template (please do not modify the structure of this file)
2. Select "Enable automatic refresh" if applicable

NOTE: The red corners give you additional information. Please do not remove the examples

3. Fill the table, beginning from the left to the right. This will display dedicated lists according to your choice. Please always use these drop down menus.
4. Select the Product Type (Custom Oligos, Double Dye Probes, Molecular Beacons, Custom siRNA, HGS Oligos, MGB Double-Dye Probes)
5. Enter your Oligonucleotides with the corresponding specifications (Green: Mandatory, Blue: Optional & White: No specification)
6. Save this Excel file on your computer
7. Upload this Excel file with your Oligonucleotides

New version 8.0

Batch file

Please note that

- Chemistry codes are:
  - ▶ ACGTA = DNA
  - ▶ (ACGLIA) = RNA
  - ▶ (ACGLIA) = 2'-O-Me RNA
  - ▶ (ACGTA) = LNA
  - ▶ A\*(G\*(U\*))A = Phosphorothioates
- Name(s) will be generated automatically if the corresponding cell is empty in your Excel file.

# HOW TO PAY

## POSTPAID SYSTEM

### One order / one invoice

You place an order of 1 or multiple oligonucleotides and you receive the invoice corresponding to this order.

## PREPAID SYSTEM Oli&GO™

### One invoice for multiple orders

You place a defined amount on your Oli&GO™ account. You receive an invoice corresponding to this amount. You can use this amount over time.



ECONOMIC

Exclusive Oli&GO™ prices.



EASY

Only one invoice for multiple oligo's orders spread over time.



CONVENIENT

Oligo orders scheduled and tracked on line.



SECURED

One administrator can give restricted access to multiple users.


The online tracking allows you to check the statements of your oligonucleotide orders at any time. >Note



# SHIPMENT GROUP


## HOW TO REDUCE MY SHIPPING FEES

With the shipment group option, all the labs from the same institution can group their shipments to benefit from free (or reduced) shipping cost.




**ECONOMIC**

Free shipping of your oligos<sup>1</sup>.



**ECOLOGIC**

Reduction of the number of parcels sent.



**FAST**

Shipping of the oligos as soon as they are ready<sup>2</sup>.

## RELATED PRODUCTS

<b>Custom genes</b>		
<b>dNTPs</b>		
dNTP Mix	1x 20 µmoles	NU-0010-10
dNTP Set	4 x 25 µmoles each	NU-0020-50
<b>Takyon™ qPCR kits</b>		
Test your free sample, visit <a href="http://www.eurogentec.com/qpcr-takyon.html">www.eurogentec.com/qpcr-takyon.html</a>		
<b>DNA purification kits</b> (100 preps)		
SMARTPure PCR Kit		SK-PCPU-100
SMARTPure Gel Kit		SK-GEPU-100
SMARTPure Plasmid Kit		SK-PLPU-100
<b>DNA extraction kit</b> (100 preps)		
SMARTEExtract DNA kit		SK-DNEX-100
<b>Agarose</b>		
Agarose - 500g		EP-0010-05
Agarose small fragment (125g)		EP-0020-10
AgaTabs - 300 tablets		EP-0030-15
<b>MW markers</b>		
SmartLadder	(200 to 10000 bp)	MW-1700-10
SmartLadder SF	(100 to 1000 bp)	MW-1800-04
<b>Electrophoresis devices</b>		
Mupid®-One	EU cable	MU-0041-
	UK cable	MU-0041+
SmartViewer for Mupid®		MU-0101
SmartIlluminator		MU-0201

With the realtime integrated counter, you keep an eye on your budget.

Scheduling your oligonucleotide orders allows reducing the number of parcels sent and decreases your shipping costs.

Functionalities	Administrator	User
	Full of privileges to control the system	Restricted access
Use 1 or more Accounts	✓	✓
Buy Oligonucleotides	✓	✓
Receive an Order Confirmation	✓	✓
Receive the Related Documentation	✓	✓
Rename the Account	✓	✗
Add/Remove User(s)*	✓	✗
Define/Update Shipping Address	✓	✗
Reload the Account(s)	✓	✗
Schedule Orders (Day/Time)	✓	✗

\*Multiple users can be defined per account

>Note  
<sup>1</sup>Minimum order amount required. <sup>2</sup>Depending on the order quantity, we can determine a delivery plan.  
 See our detailed shipping conditions on <http://www.eurogentec.com/shipping-conditions.html>



1. RNA interference
  - RNA interference is licensed under European Patents T144623, I214945 and foreign equivalents from Alnylam Pharmaceuticals, Inc., Cambridge, USA and is provided only for use in academic and commercial research (excluding the evaluation or characterization of this product as the potential basis for a siRNA-based drug) and not for any other commercial purposes. Information about licenses for commercial use (including discovery and development of siRNA-based drugs) is available from Alnylam Pharmaceuticals, Inc., 300 Third Street, Cambridge, MA 02142, USA.
2. Elitech Dyes and Quenchers
  - Use of this product (including but not limited to Yakima Yellow® and Eclipse™ Dark Quencher) is covered by patents owned by ELITechGroup, Inc., including but not limited to the following US patents and corresponding patent claims outside the US: US 6,699,975, US 6,790,945, US 7,662,942, US 6,653,473, US 6,972,339, 7,541,454, 7,671,218, 7,767,834, 8,163,910, 6,727,356. This product, including both Research Grade Product and IVD Grade Product, is sold under a licensing arrangement with ELITechGroup, Inc. for the purchaser's own internal research. Commercial use of this product requires a separate license from Kaneka Eurogentec S.A. Further information on purchasing licenses may be obtained by contacting [licensing@eurogentec.com](mailto:licensing@eurogentec.com).
3. Alexa Fluor®, BODIPY®, Cascade Blue®, Marina Blue®, Oregon Green®, Pacific Blue®, QSY®, Rhodamine Red™ Texas Red® dyes
  - This product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is conditioned on the buyer using the purchased product solely in research conducted by the buyer, excluding contract research or any fee for service research, and the buyer must not (1) use this product or its components for (a) diagnostic, therapeutic or prophylactic purposes; (b) testing, analysis or screening services, or information in return for compensation on a per-test basis; or (c) manufacturing or quality assurance or quality control, and/or (2) sell or transfer this product or its components for resale, whether or not resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5781 Van Allen Way, Carlsbad, CA 92008 USA or [outlicensing@thermofisher.com](mailto:outlicensing@thermofisher.com).
4. Black Hole Quencher<sup>†</sup>
  - Black Hole Quencher® and BHQ® are registered trademarks of Biosearch Technologies, Inc., Novato, California, U.S.A. ('BTI') Patents are currently pending for the BHQ technology and such BHQ technology is licensed by the manufacturer pursuant to an agreement with BTI and these products are sold exclusively for research and development use only. They may not be used for human or veterinary in vitro or clinical diagnostic purposes and they may not be re-sold, distributed or re-packaged. For information on licensing programs to permit use for human or veterinary in vitro or clinical diagnostic purposes, please contact Biosearch at [licensing@biosearchtech.com](mailto:licensing@biosearchtech.com).
5. C-5 Propyne and G-clamp modified Oligonucleotides
  - C-5 Propyne and G-clamp modified oligonucleotides are covered by patents or patents pending owned by Isis Pharmaceuticals, Inc. ("Isis"). Purchase of this product includes a limited license to use this product solely for internal research. This license specifically excludes (and you have no right to use this product for): (a) therapeutic or diagnostic applications (including products or services that incorporate this product), (b) any in vivo toxicity/safety study in support of an investigational new drug application (or foreign counterpart), (c) resale (including sale of any products or services that incorporate this product) or (d) gene functionalization activities (including products or services that incorporate data derived from gene functionalization activities) if such activities have commercial application, any and all of which require a separate license from Isis. Neither this product nor any product created through its use may be used in human clinical trials. In the event you have a separate agreement with Isis regarding this product which explicitly states that the foregoing is not applicable to you, your use of this product will be governed by the terms of such agreement. In no event does the limited license included with the purchase of this product expand or alter the scope of the license granted pursuant to such agreement.
6. PNA
  - All custom PNA oligomers sold by Kaneka Eurogentec are manufactured by Panagene under an exclusive license and to be used for internal research only. And such custom PNA oligomers are not to be resold unless properly licensed to do so by separate license. All custom PNA oligomers are made and sold pursuant to license under one or more of US Patents Nos. US 5,773,571, 6,133,444, US 6,172,226, US 6,395,474, US 6,414,112, US 6,613,873, US 6,710,163 and US 6,713,602, or corresponding patent claims outside the US.
7. Hilyte™ Fluor Dyes
  - Use of this product is covered by one or more of the following US patents owned by AnaSpec, Inc.: USP 7,465,810, USP 7,754,893, USP 7,820,783, USP 8,258,292, and divisionals, continuations, continuations-in-part, reissues, substitutes, and extensions thereof. The purchase of this product conveys to the buyer a limited, non-exclusive, nontransferable right (without the right to sell, repackage, or further sublicense) under these patent rights to use only the amount of product purchased for the purchaser's own internal research. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of this product does not include nor carry any right or license to use, develop, or otherwise exploit this product commercially, and no rights are conveyed to the buyer to use the product or components of the product for any other purposes, including without limitation, provision of services to a third party, generation of commercial databases, or clinical diagnostics or therapeutics. This product is for research use only. For information on purchasing a license to this product for purposes other than research, contact [licensing@eurogentec.com](mailto:licensing@eurogentec.com).
8. QXL<sup>®</sup> Quenchers
  - Use of this product is covered by one or more of the following US patents owned by AnaSpec, Inc.: USP 7,910,753, USP 8,093,411, and divisionals, continuations, continuations-in-part, reissues, substitutes, and extensions thereof. The purchase of this product conveys to the buyer a limited, non-exclusive, nontransferable right (without the right to sell, repackage, or further sublicense) under these patent rights to use only the amount of product purchased for the purchaser's own internal research. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of this product does not include nor carry any right or license to use, develop, or otherwise exploit this product commercially, and no rights are conveyed to the buyer to use the product or components of the product for any other purposes, including without limitation, provision of services to a third party, generation of commercial databases, or clinical diagnostics or therapeutics. This product is for research use only. For information on purchasing a license to this product for purposes other than research, contact [licensing@eurogentec.com](mailto:licensing@eurogentec.com).
9. Minor Groove Binder and Eclipse™ Dark Quencher
  - Use of this product is covered by patents owned by ELITechGroup, Inc., including but not limited to the following US patents and corresponding patent claims outside the US: US 6,492,346, US 6,884,584, US 11,784,539, US 7,556,923, US 12,172,999, US-2009-0048427-A1, US 6,084,102, US 6,426,408, US 6,727,356, US 6,790,945, US 6,699,975, US 6,486,308 and US 7,794,945. This product, including both Research Grade Product and IVD Grade Product, is sold under a licensing arrangement with ELITechGroup, Inc. for use by purchaser only in human in vitro diagnostics applications to provide information for individual patient management solely using 5' nuclease (Tagman-MGB) technology. Commercial use of this product requires a separate license from Kaneka Eurogentec S.A. Further information on purchasing licenses may be obtained by contacting [licensing@eurogentec.com](mailto:licensing@eurogentec.com).
10. MOE
  - The use of this product is covered by patents or patents pending owned by Ionis Pharmaceuticals, Inc. ("Ionis"). Purchase of this product includes a limited, non-transferable right to use the purchased amount of the product for the purchaser's own internal research and does not include any license to use this product for other purpose like (for instance and not exhaustive): (a) therapeutic or diagnostic applications, (b) any in vivo toxicity/safety study in support of an investigational new drug application (or foreign counterpart), (c) resale (including sale of any products or services that incorporate this product or oligonucleotides) or (d) gene functionalization activities (including products or services that incorporate data derived from gene functionalization activities). Any and all of such application of the oligonucleotides may require a separate license from Ionis.

## EUROGENTEC TRADEMARKS AND LABELS

Alexa Fluor® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)

BHQ-0®, BHQ-1®, BHQ-2®, BHQ-3® are registered trademarks of Biosearch Technologies, Inc.

Black Hole Quencher® is a registered trademark of Biosearch Technologies, Inc.

CAL Fluor™ is trademark of Biosearch Technologies, Inc.

BODIPY® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)

Cascade Blue® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)

Cy® is a registered trademark of GE Healthcare Bio-Sciences Corp.

Cy®5 Direct is a registered trademark of GE Healthcare Bio-Sciences Corp.

Dragonfly Orange™ is a trademark of Kaneka Eurogentec S.A.

Eclipse® is a registered trademark of Epoch Biosciences, Inc.

Excel™ is a trademark of Microsoft Corp.

HiLyte™ is a trademark of Anaspec, Inc.

LNA® is a registered trademark of Exiqon A/S.

Marina Blue® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)

MGB® is a registered trademark of Epoch Biosciences, Inc.

MGB Eclipse® is a registered trademark of Epoch Biosciences, Inc.

Oli&GO™ is a trademark of Kaneka Eurogentec S.A.

Oregon Green® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)

Pacific Blue™ is a trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)

Plexor™ is a trademark of Promega Corporation

Quasar™ is trademark of Biosearch Technologies, Inc.

QXL® is a registered trademark of Anaspec, Inc.

Rhodamine Red™ and Rhodamine Green™ are trademarks of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)

ROX™ is a trademark of Applera Corp.

RP-Cartridge-Gold™ is a trademark of Kaneka Eurogentec S.A.

Taqman® is a registered trademark of Applied Biosystems, Inc.

Texas Red® is a registered trademark of Molecular Probes, Inc. (a wholly owned subsidiary of Invitrogen Corp.)

Track™ is a trademark of Kaneka Eurogentec S.A.

Unique™ is a trademark of Kaneka Eurogentec S.A.

UltraPureGold™ is a trademark of Kaneka Eurogentec S.A.

Vic® is a registered trademark of Applera Corp.

Yakima Yellow® is a registered trademark of Epoch Biosciences, Inc.

www.eurogentec.com

# Oligos

**Contact:**

For more information,  
please contact your Sales  
representative or contact us at:  
[oligo@eurogentec.com](mailto:oligo@eurogentec.com)

**Kaneka Eurogentec S.A.**

rue Bois Saint-Jean, 5  
4102 Seraing  
Belgium

Eurogentec 2019

Follow us on



## Life Science Research



Oligos




PCR/qPCR



Genes/DNA



Antibodies



Peptides



Proteins